


PROCEDURE



PART 1: EXTRACTION OF CAFFEINE:

- Measure out 100 ml of soda pop.
- Crush the sample with mortar and pestle and dissolve in 100 ml of water.
- Add 2 g of sodium carbonate and dissolve.
- Transfer the solution into a separatory funnel.
- In a hood, add 10 ml of methylene chloride.
- Using silicone put the lid and rotates.
- Mix until no more gas is evolved.
- Decant the organic layer into a 50 ml Erlenmeyer flask.
- Add 8 ml of methylene chloride to the aqueous layer, mix and add the organic layer to the first extract. Repeat again this extraction.
- Add 1 g of anhydrous Na₂SO₄.
- Wait 10 min.
- Weight a 50 ml erlenmeyer flask.



Part 2: PURIFICATION OF CAFFEINE:

- Transfer the extract into the preweighed flask
- Set up for simple distillation and drive off most of the solvent.
- Purify the caffeine by sublimation.



Part 3: TLC PROCEDURES:

- Cool the flask and reweight it to determine the amount of caffeine.
- Use a 250 ml beaker lined with a piece of filter paper covered with a watch glass.
- Prepare the developing chamber by obtaining a small amount of the development solvent
- Fill the chamber to a depth of about ¼ inch.
- Draw a line with a pencil on a TLC plate and mark off two intervals equally spaced on the line.
- In the hood add 5 ml of methylene chloride in the flask and dissolve the crystals.
- Spot the solution on the plate and dip a capillary tube in the solution. Allow, just a small part of the liquid on the plate.
- When the plate has been spotted, it's ready to be going on, control the level of the solvent must be checked by applying stains on the plate because they will dissolve on it inside the chamber
- When the solvent has risen to about 1 centimeter on the end of the slab, mark the position of the spots with a pencil (surrounding the spots).
- When the plate is dry, check the stains with a UV lamp and observe all the points formed, and surround them with a pencil



Part 4: TLC PROCEDURES:

- Make the melting point of the obtained caffeine and compare the value obtained with the original one.